Mitra® Device (IVD)

These instructions apply to the following Mitra® Device configurations:

- Mitra Device Alone
- Mitra Device in Specimen Bag with Desiccant


IMPORTANT – Read the entire instructions before use!

Intended Use:

A single-use, non-sterile device used as a specimen collector, and for the storage and transport of blood and other biological fluids for analytical and diagnostic analyses.

Device Description:

The intended user is a healthcare professional or a layperson. The blood sample is typically taken through a prick in the fingertip using a lancet. The sampler tip (contained within the Mitra Device) is gently applied to the drop of blood on the fingertip until the entire sampler tip turns red indicating a volumetric sample (10, 20, or 30 μ L depending on the device being used) has been collected. This process is repeated for each sampler tip in the device. Once all sampler tips (one (1), two (2), three (3), or four (4) depending on the device being used) have been filled, the device outer housing is clicked closed to secure the sample(s).

Materials Provided:

- Mitra Device
- Specimen Bag with Desiccant
- Sample ID Barcode

Materials Required, Not Provided by Trajan:

- Disposable Lancet**
- Gauze
- Bandage
- Specimen Bag with Desiccant*

* Part #s containing an “A” represent the “alone” configuration and do not contain a specimen bag (e.g. VM-102ANLR). Part #s containing a “S” represent the “specimen bag” configuration and do contain a bag (e.g. VM-102SNLR).


** Use a contact-activated lancet with a 1.5 mm wide and 2 mm deep cut to obtain sufficient blood without squeezing for adult populations. For vulnerable populations use a contact-activated lancet with a 1 mm wide and 2 mm deep cut.

Product Specifications:

- Sample Type: Dried Whole Blood
- Formats/Configurations: Alone or In Specimen Bag with Desiccant
- Sample Volume: 10, 20, or 30 μ L
- Substrate: Hydrophilic porous polymer
- Volumetric Precision (%RSD): < 5%
- Number of Samples Collected: Up to 4 (depending on Device used)
- Shelf Life: 18 months
- Device Storage: Up to 30°C. For a shorter time (e.g. during transport) up to 45°C.
- Sample Storage: The guidelines for sample storage are analyte dependent and will need to be determined by device end-user.

Troubleshooting:

- Not enough blood was produced:
 - Warm hands thoroughly (e.g. run hands under hot water, rub hands together, walk briskly for 5 minutes).
 - Shake and massage fingers to activate blood flow.
 - Prick finger with a new lancet.
- Sampler tip is not completely red:
 - The tip was removed from the blood too soon. Reapply the tip to the blood drop until no white remains visible on sampler tip.
 - Blood flow stopped. See suggestions above.
 - Blood was dropped onto the device tip from above. Always point sampler tip towards floor when touching it to the blood drop.

Warnings and Precautions:

- Single use only for a single individual.
- Do not use after expiration date.
- For external use only.
- Devices should be transported/mailed to the analytical laboratory, and appropriate documentation maintained according to local regulations and the analytical laboratory procedures and policies.
- Do not use if device packaging has been opened or damaged.
- Labs must validate use of product for their specific assay.
- Observe universal biological risk precautions.
- All used materials with blood residues must be handled and disposed of safely in accordance with local regulations.
- Any serious incident that has occurred in relation to the medical device should be reported to Trajan Scientific and Medical and the competent authority of the Member State in which the user and/or patient is established.

WATCH – Instructional videos at www.neoteryx.com/collect or scan the QR code to access.

Instructions for Sampling:

1.OPTIONAL: If your device came in a silver specimen bag, tear open the specimen bag using the pre-slit notch. Remove the device from the specimen bag. Do not remove the desiccant from specimen bag. If your device did not come in a bag, skip to step 2.

2.Open device by pulling tabs on device apart from each other until the two flaps are folded downwards.

NOTE: Do NOT remove sampler tips from samplers. Do not remove samplers from device.

3.Touch sampler tip to surface of blood sample.


Watch sampler tip turn FULLY red as it absorbs the blood, count 2 seconds, and remove. It is OK to apply the sampler tip to the blood drop many times to fill until no white is visible.

IMPORTANT!

- Do NOT fully submerge the sampler tip in blood.
- Sampler tip should always point downward towards floor as illustrated.
- Do NOT drip blood onto the sampler tip.

4.Repeat step 3 with remaining samplers in device. Device will contain one (1), two (2), three (3), or four (4) samplers depending on device format received.


5.Check the sampler tips. Upon successful sampling, the entirety of each sampler tip should be colored red with no white visible.


Ensure sampler tips are filled correctly.

Under-sampling occurs when:

1. The sampler tip is removed from the blood too soon. Touch tip to blood until no white remains.
2. Blood flow stops. If this happens, repeat step 3 with same tip until it turns fully red.

Over-sampled

Under-sampled


Correctly sampled

Over-sampling occurs when blood is dripped onto the sampler tip from above.

6.Lift both flaps of the device up to meet at the top and press tabs together until a click is heard.

7.IF YOU RECEIVED A SPECIMEN BAG:
Samplers secured in device should be immediately placed in provided specimen bag containing desiccant and SEALED tightly using the zipper-top closure.

**For further assistance,
please do not hesitate to contact us:
neo.support@trajanscimed.com**

Disclaimer:

Mitra® devices are CE-IVD (IVDR) devices intended as specimen collectors and for the storage and transport of dried blood. They are available as registered IVD Devices in the European Union and United Kingdom, Australia, Brazil, China, and Canada, as well as multiple Health Ministries worldwide. In the USA, Mitra devices are supplied as a research use only (RUO) product to assist in method development, other research-related and non-diagnostic activities. Mitra has not been validated for use with any diagnostic testing.

Specifications are subject to change without notice.

Neoteryx®, Mitra® and VAMS® are registered trademarks owned by Trajan Scientific Australia Pty Ltd.

 Trajan Scientific Australia Pty Ltd (Sponsor) 7 Argent Place, Ringwood, VIC 3134, AU - Tel: +61 (0) 3 9874 8577		
EC	REP	Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta
 www.neoteryx.com/collect		

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Neoteryx®, Mitra® e VAMS® são marcas registadas propriedade da Trajan Scientific Australia Pty Ltd.