Mitra® Device Collection Kit (Procedure Pack)

These instructions apply to the following Mitra® Device configurations:

- Microsample Collection Kit


IMPORTANT – Read the entire instructions before use!

Intended Use:

The Microsample Collection Kit is a single-use non-sterile procedure pack intended for the collection, storage and transport of capillary whole blood specimens of specified volume using the contents of the collection kit.

Device Description:

The intended user is a healthcare professional or a layperson. The blood sample is typically taken through a prick in the fingertip using a lancet. The VAMS® sampler tip (contained within the Mitra Device) is gently applied to the drop of blood on the fingertip until the entire sampler tip turns red indicating a volumetric sample (10, 20, or 30 µL depending on the device being used) has been collected. This process is repeated for each sampler tip in the device. Once all sampler tips (one (1), two (2), three (3), or four (4) depending on the device being used) have been filled, the device outer housing is clicked closed to secure the sample(s).

Materials Provided:


- Mitra Device
- Specimen Bag with Desiccant
- Sample ID Barcode
- Lancet (x2 or x4)
- Gauze (x1 or x2)
- Bandage (x2 or x4)
- Plastic Shipping Envelope

Product Specifications:

- Sample Type: Dried Whole Blood
- Formats/Configurations: Collection Kit
- Sample Volume: 10, 20, or 30 µL
- Substrate: Hydrophilic porous polymer
- Volumetric Precision (%RSD): < 5%
- Number of Samples Collected: Up to 4 (depending on Device used)
- Shelf Life: 18 months
- Device Storage: Up to 30°C. For a shorter time (e.g. during transport) up to 45°C.
- Sample Storage: The guidelines for sample storage are analyte dependent and will need to be determined by device end-user.

Troubleshooting:

- Not enough blood was produced:
 - Warm hands thoroughly (e.g. run hands under hot water, rub hands together, walk briskly for 5 minutes).
 - Shake and massage fingers to activate blood flow.
 - Prick finger with a new lancet.
- Sampler tip is not completely red:
 - The tip was removed from the blood too soon. Reapply the tip to the blood drop until no white remains visible on sampler tip.
 - Blood flow stopped. See suggestions above.
 - Blood was dropped onto the device tip from above. Always point sampler tip towards floor when touching it to the blood drop.

Warnings and Precautions:

- Single use only for a single individual.
- Do not use after expiration date.
- For external use only.
- Devices should be transported/mailed to the analytical laboratory, and appropriate documentation maintained according to local regulations and the analytical laboratory procedures and policies.
- Do not use if device packaging has been opened or damaged.
- Labs must validate use of product for their specific assay.
- Observe universal biological risk precautions.
- All used materials with blood residues must be handled and disposed of safely in accordance with local regulations.
- Any serious incident that has occurred in relation to the medical device should be reported to Trajan Scientific and Medical and the competent authority of the Member State in which the user and/or patient is established.

WATCH – Instructional videos at www.neoteryx.com/collect or scan the QR code to access.

Instructions for Sampling:

1. Unpack and prepare your kit contents:

a. Twist off lancet caps.

b. Tear open silver specimen bag.

c. Remove device from specimen bag.

IMPORTANT – Leave the desiccant in the specimen bag!

d. Open the device by pulling tabs on device apart from each other until the two flaps are folded downwards. **NOTE:** Do NOT remove sampler tips from sampler bodies. Do NOT remove sampler bodies from device.

e. Unwrap gauze.

f. Unpeel back of bandages.
2. **Wash hands** with warm, soapy water.




3. **Warm your hands** by rubbing them together. You may also warm your hands by walking, jumping, or any other method you prefer. **Warming your hands is IMPORTANT for good blood flow.**

4. **Select the finger you will use** to collect the blood sample. See illustration for recommended locations off-center of finger midline.

5. **Prick your finger** by laying hand on a hard surface. Position the lancet on the chosen location. Apply firm pressure to end of lancet until a click is heard.

6. **Wipe away first blood drop with the gauze.** A blood drop may take up to 20 seconds to form. If needed, gently massage pricked finger upward until drop forms.

7. **Collect sample by touching one of the sampler tips to the blood drop.** First, watch the tip turn fully red as it absorbs the blood. Next, count 2 seconds. Last, slowly remove the sampler tip from blood. Repeat this process with the remaining sampler tips.


You may apply sampler tips to the blood drop multiple times to fill. Always point sampler tip towards floor when touching it to blood drop on your finger as shown above. Do NOT drip blood onto the sampler tip.


Ensure sampler tips are filled correctly.

- Under-sampling occurs when:
1. The sampler tip is removed from the blood too soon. Touch tip to blood until no white remains.

2. Blood flow stops. If this happens, repeat steps 4-7 with same tip until it turns fully red.

Over-sampled

Under-sampled

Correctly sampled

Over-sampling occurs when blood is dripped onto the sampler tip from above. Always touch sampler tip to blood drop as shown in step 7.



8. **Set device down** on hard surface and apply bandage to your finger.

If you experience excessive bleeding after collection, apply pressure and contact your healthcare provider immediately.
9. **Close device** by lifting the flaps to meet at the top. Press together until a click is heard.

10. **Insert device** into the specimen bag and seal shut. **Ensure the desiccant is still in the bag.**

11. **Insert specimen bag** into plastic shipping envelope. Ensure a pre-addressed label is affixed to mailing envelope. Send to lab for processing as instructed.

12. **Dispose of the lancet** by enclosing it in the plastic pouch that contained the lancets, seal it and trash it in a bin. Do not dispose the lancet in a recycling bin.

For further assistance, please do not hesitate to contact us:
neo.support@trajanscimed.com

Disclaimer:
Mitra® devices are CE-IVD (IVDR) devices intended as specimen collectors and for the storage and transport of dried blood. They are available as registered IVD Devices in the European Union and United Kingdom, Australia, Brazil, China, and Canada, as well as multiple Health Ministries worldwide. In the USA, Mitra devices are supplied as a research use only (RUO) product to assist in method development, other research-related and non-diagnostic activities. Mitra has not been validated for use with any diagnostic testing.

Specifications are subject to change without notice.

Neoteryx®, Mitra® and VAMS® are registered trademarks owned by Trajan Scientific Australia Pty Ltd.

		
Trajan Scientific Australia Pty Ltd (Sponsor) 7 Argent Place, Ringwood, VIC 3134, AU - Tel: +61 (0) 3 9874 8577		
EC REP		Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta
		www.neoteryx.com/collect

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Neoteryx®, Mitra® e VAMS® são marcas registadas propriedade da Trajan Scientific Australia Pty Ltd.