

Harpera Microbiopsy Tool (IUO): Safety Statement

Harpera Device Description and Intended Use

Trajan Harpera (IUO) is a medical device currently under development for the collection of a microbiopsy from skin, for use in clinical studies. In the hands of a trained healthcare professional, this minimally invasive device utilizes a novel technique to collect skin biopsy specimens, including from sensitive areas around the face. This single-use device is supplied individually wrapped and has been exposed to a gamma irradiation process similar to processes utilized for the sterilization of other surgical instruments. Harpera is not associated with any specific diagnostic test or testing scheme and is not intended for the analyses of a specific disease states.

Data From Clinical Studies

Several clinical studies have been conducted by Trajan collaborators using Harpera prototype devices. No adverse events were reported to Trajan and no field safety corrective actions (FSCA) were requested during the course of these studies.

Several studies investigated subjects' pain response during the biopsy process utilizing Harpera. A pain score was recorded and the results indicate that study subjects did not experience significant pain during the procedure [*Lin et al. (2013)*, *Churiso et al. (2020)*, *Yamada et al. (2020)*].

Currently, a Human Factor Engineering study is being planned to evaluate the safety parameters for the use of Harpera. Included will be a survey to determine ease-of-use of the device as well as comprehension of the instructions for use (IFU).

Regulatory Status

Harpera is a Development Stage Device and is supplied as "Investigational Use Only" (IUO) for clinical studies, with appropriate labeling. Validation of the design, performance and safety of the device has not been completed.

Stability of a Harpera skin specimen has been investigated through genomic techniques (i.e. detection of common housekeeping RNA, human DNA). However, validation of specimen stability for specific biomarkers is beyond the intended purpose of the device.

Trajan plans to initiate the regulatory filings in several jurisdictions in 2024.

Anne Marie Whalen, PhD
Global Head, Regulatory Affairs.

*Electronically signed
by: Anne Marie Whalen
Reason: Regulatory
representative signs to
confirm document
complies with
applicable regulatory
requirements.
Date: Jan 16, 2024
13:30 EST*

References

References of the Harpera prototype in human studies:

- Lin, L. L., et al. (2013). "Microbiopsy engineered for minimally invasive and suture-free sub-millimetre skin sampling." [F1000Res 2: 120.](#)
- Banan, P., et al. (2013). "Effects of ex vivo skin microbiopsy on histopathologic diagnosis in melanocytic skin lesions." [JAMA Dermatol 149\(9\): 1107-1109.](#)
- Prow, T. W., et al. (2013). "The opportunity for microbiopsies for skin cancer." [Future Oncol 9\(9\): 1241-1243](#)
- McClenahan, P., et al. (2014). "BRAFV600E mutation status of involuting and stable nevi in dabrafenib therapy with or without trametinib." [JAMA Dermatol 150\(10\): 1079-1082](#)
- Tan, J.-M., et al. (2015). "BRAFWild-Type Melanoma in Situ Arising In a BRAFV600E Mutant Dysplastic Nevus." [JAMA Dermatology 151\(4\).](#)
- Dermatology Research Centre, School of Medicine, The University of Queensland (March 2016), "[Natural history and properties of naevi in advanced melanoma patients receiving treatment](#)", ACTRN 12616000272493 (ANZCTR)
- Tom, L. N., et al. (2016). "Skin microbiopsy for HPV DNA detection in cutaneous warts." [J Eur Acad Dermatol Venereol 30\(12\): e216-e217.](#)
- Sobarun, P., et al. (2017). "Microbiopsy Biomarker Profiling in a Superficial Melanoma Resembling a Pigmented Basal Cell Carcinoma." [JAMA Dermatol 153\(4\): 334-336.](#)
- Kirstein, O. D., et al. (2017). "Minimally invasive microbiopsies: a novel sampling method for identifying asymptomatic, potentially infectious carriers of Leishmania donovani." [Int J Parasitol 47\(10-11\): 609-616.](#)
- Lin, L. L. and T. W. Prow (2017). "Novel microdevices for controlled blood and skin extraction, NHMRC." [Impact 2017\(6\): 58-60.](#)
- Lei, B. U. W., et al. (2019). "Absorbent Microbiopsy Sampling and RNA Extraction for Minimally Invasive, Simultaneous Blood and Skin Analysis." [JoVE\(144\): e58614.](#)
- Churiso, G., et al. (2020). "Minimally Invasive Microbiopsies as an Improved Sampling Method for the Diagnosis of Cutaneous Leishmaniasis." [Open Forum Infect Dis 7\(9\): ofaa364.](#)
- Yamada, M., et al. (2020). "A minimally invasive clinical model to test sunscreen toxicity based on oxidative stress levels using microbiopsy and confocal microscopy - A Proof of concept study." [Int J Cosmet Sci.](#)
- Yamada, M., et al. (2020). "Microbiopsy-based minimally invasive skin sampling for molecular analysis is acceptable to Epidermolysis Bullosa Simplex patients where conventional diagnostic biopsy was refused." [Skin Res Technol.](#)
- Michele Fimiani, P. R., Elisa Cinotti (2020). Technology in Practical Dermatology, [Springer Cham.](#)
- Cloots, K., et al. (2021). "Assessing L. donovani Skin Parasite Load: A Proof of Concept Study of a Microbiopsy Device in an Indian Setting." [Front Cell Infect Microbiol 11: 645121.](#)
- Owen, S. I., et al. (2021). "Evaluation of qPCR on blood and skin microbiopsies, peripheral blood buffy coat smear, and urine antigen ELISA for diagnosis and test of cure for visceral leishmaniasis in HIV-coinfected patients in India: a prospective cohort study." [BMJ Open 11\(4\): e042519.](#)
- Haderl E, Mosca M, Hong J, Brownstone N, Liao W, Bhutani T. Innovations in translational research in dermatology: minimally invasive methods for biosample acquisition. [Dermatol Online J. 2021 Oct 15;27\(10\).](#)
- Preis, S., et al. (2022). "Munich atopy prediction study (MAPS): protocol for a prospective birth cohort addressing clinical and molecular risk factors for atopic dermatitis in early childhood." [BMJ Open 12\(9\): e059256.](#)
- Jain, M., et al. (2022). "Minimally invasive microbiopsy for genetic profiling of melanocytic lesions: A case series." [J Am Acad Dermatol 87\(4\): 903-904.](#)
- Liverpool School of Tropical Medicine (April 2023). "[Using BCG Vaccine to Understand Tuberculosis Infection](#)". [NCT05820594 \(ClinicalTrials.gov\)](#)
- [Product Specification Sheet for Harpera](#)

Safety statement Harpera IUO_rev2

Final Audit Report

2024-01-16

Created:	2024-01-15 (Central European Time)
By:	Florian Lapierre (flapierre@trajanscimed.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAf9O5BWMzoKIljCnpWxXICDaFdH4-gsKY

"Safety statement Harpera IUO_rev2" History

-  Document created by Florian Lapierre (flapierre@trajanscimed.com)
2024-01-15 - 17:19:40 GMT+1
-  Document emailed to Anne Marie Whalen (AWhalen@trajanscimed.com) for signature
2024-01-15 - 17:20:15 GMT+1
-  Email viewed by Anne Marie Whalen (AWhalen@trajanscimed.com)
2024-01-15 - 17:20:21 GMT+1
-  Anne Marie Whalen (AWhalen@trajanscimed.com) authenticated with Adobe Acrobat Sign.
Challenge: The user opened the agreement.
2024-01-16 - 19:13:47 GMT+1
-  Anne Marie Whalen (AWhalen@trajanscimed.com) authenticated with Adobe Acrobat Sign.
Challenge: The user completed the signing ceremony.
2024-01-16 - 19:30:41 GMT+1
-  Document e-signed by Anne Marie Whalen (AWhalen@trajanscimed.com)
Signing reason: Regulatory representative signs to confirm document complies with applicable regulatory requirements.
Signature Date: 2024-01-16 - 19:30:43 GMT+1 - Time Source: server
-  Agreement completed.
2024-01-16 - 19:30:43 GMT+1



Adobe Acrobat Sign