TRAJAN

Trajan Scientific and Medical

Overview

We present four studies associated with blood microsampling using a Dried Blood Spot (DBS) workflow on paper and polymer substrates:

- 1. Targeted metabolomics analysis of 2.74 μL blood biological and technical replicates using GC-MS/MS.
- 2.Comparison of sub-punch and accurate volume DBS for amino acids and acyl-carnitines across of range of hematocrits using FIA-MS/MS.
- 3.Development of a synthetic polymer DBS substrate for blood derived FAME analysis by GC-MS.
- 4.Extraction efficiency of DNA from different DBS substrates a comparison of paper and synthetic polymer.



Figure 1. The hemaPEN[®] has been designed to collect replicate, volumetrically accurate (2.74 µL) and precise volumes of blood. Shown here collecting a capillary blood sample, it could be used with any blood source. An embedded desiccant is designed to dry the collected sample inside the contained device. The device does not contain a lancet.

Introduction

- Blood microsampling, such as dry blood spotting (DBS) is safe, convenient, "patient centric" and ideal for expansive population health initiatives as well as research-focused "omics" studies.
- DBS adoption has been limited by the perceived lack of volumetric accuracy and the analytical bias associated with variable blood hematocrit (HCT) and analyte extraction from the paper substrate.
- This work introduces a novel device for collecting very small blood volumes (<3 µL) and their subsequent metabolite and DNA analysis of the dried blood extracts from paper substrates as well as a panel of novel porous polymer substrates. The substrates are encapsulated in a blood collection device designed to collect an accurate volume and stored in a secure desiccated environment.

Incorporating novel synthetic DBS substrates into a blood microsampling device for multi-omics analyses.

Kyle Bachus¹, Ricardo Neto⁶, Lada Staskova^{2, 3}, Konstantinos Kouremenos^{1,7}, Jeff Craig^{3, 4}, Robert Shellie⁵, Dario Arrua⁶, Emily F Hilder⁶, Andrew Gooley^{1,7}, Wei Boon Hon¹ 1) Trajan Scientific and Medical, Ringwood, Australia, 2) RMIT University, Melbourne, Australia, 3) Centre for Molecular and Medical Research, School of Medicine, Faculty of Health Deakin University, Geelong, Australia, 4) Murdoch Children's Research Institute, The Royal Children's Hospital, Parkville, Australia, 5) Centre for Advanced Sensory Science (CASS), School of Exercise and Nutrition Sciences, Deakin University, Melbourne, Australia, 6) Future Industries Institute, University of South Australia, Mawson Lakes Campus, Adelaide, Australia, 7) Bio21 Institute, University of Melbourne, Australia.

Methods

Study 1

Three biological replicates were collected within the same time point (within 5 min) from a healthy donor by fingerstick capillary blood. Each sample collects $4 \times 2.74 \mu$ L technical replicates of whole blood onto pre-punched grade PerkinElmer 226^{TM} filter paper. Samples were extracted and analyzed as TMS derivatives by GC-MS/MS. Sample preparation - each technical replicate (12 samples in total) was extracted into 100 μ L of chilled methanol (IS at 5 μ M); dried and converted to TMS-derivatives. A 1 μ L sample was injected into a Shimadzu GCMS-TQ8040 (utilizing Shimadzu Smart Metabolites Database). Data was processed using Shimadzu's Lab Solutions software, and the matrix output was then imported into Metabo Analyst for statistical analysis. Log transformation and median normalization (with no scaling) was applied as the data pretreatment step in order to reduce the differences between possible metabolite concentrations and influence of other factors such as measurement noise.

Study 2

Preparation of blood - five different HCT levels were prepared from a single donor. Eight hemaPEN were used to collect and dry samples at each HCT for the accurate volume workflow (n = 32 pre-punch WhatmanTM 903 specimen filter paper). DBS cards (WhatmanTM 903 specimen filter paper – 4 circles each containing 75 µL of the prepared blood) were used for the sub-punch workflow where 4 × 3.2 mm sub-punches were taken from each of the blood spots (n = 32 subpunch DBS). All DBS were extracted with 150 µL of 100% MeOH containing the appropriate deuterated standard; shaken for 60 min and dried under an air stream at 65°C then reconstituted in 200 µL of 70% MeCN. 10 µL samples were flow injected using the Acquity UPLC (SM-FTN) into a Waters Acquity TQ MS detector. Twenty of the extracted metabolites (amino acids and acylcarnitines) were analysed and reported in µmol/L based on the peak ratio of the extracted analyte and corresponding deuterated IS. Data analytics on the µmol/L peak ratios was performed using SIMCA-14 (Sartorius Stedim Data Analytics).

Study 3

10 μ L and 3 μ L blood samples were applied to both 226 paper substrate (6 mm diameter) as well as the synthetic polymer mPPM inserts. Samples were dried and added to 2 mL of 1% (v/v) H²SO⁴ in anhydrous methanol. Transmethylation reaction was performed at 70°C for 3 h. FAMEs were extracted into 600 μ L of n-heptane and analysed by GC-MS (Thermo Scientific FOCUS GC with DSQ II MS). Whole blood was also transmethylated as described above and used as the reference for all DBS workflows.

Study 4

Five different paper substrates and the mPPM insert were used to explore the extraction performance using two different commercial kits Qiagen QIAamp DNA and Zymo Quick DNA mini. Four x 3 mm sub-punch samples were collected from each card and a 12 µL blood sample was collected onto the mPPM insert. All samples were incubated overnight in PBS and Proteinase K solution prior to extraction. DNA concentration was measured by Qubit Fluorometric Quantitation and high molecular weight quality determined by agarose gel electrophoresis.







R2X[2] = 0.232

R2x[1] = 0.444

Ellipse: Hotelling's T2 (95%)

Figure 3. PCA Score Plot of all 20 amino acids & acyl-carnitines peak ratios for accurate volume pre-punched DBS (Grey) and sub-punch DBS (Orange) (n=32 for each of the 5 HCT and each workflow).

PC1 - Observed variation attributed to analytes influenced by the sub-punch workflow.

PC2 - Observed variation attributed to analytes influenced by the HCT of the accurate volume samples. Inset shows the same PCA score plot highlighting the volumetric bias at low and high HCTs for sub-punch workflows compared to the accurate volume hemaPEN[®] sampling.

Figure 4. PCA Score Plot of all acyl-carnitine peak ratios for accurate volume pre-punched DBS (H-series) and sub-punch DBS (G-series) workflows for five different HCTs (n=32 for each HCT and each workflow PC1 accounts for 44.2% of the observed variation which can be attributed to the HCT regardless of the DBS workflow. PC2 accounts for 23.2% of the observed variation which can be attributed to the different DBS workflows – both the sub-punch and accurate volume are well separated at 63%, 55% and 25% HCTs.

Study 3



Figure 5. hemaPEN cartridge with mPPM inserts (3.6 mm ID x 3.5 mm depth) and 2.74 μ L capillary blood collection.



different DBS substrates.

Conclusion:

Here we have characterized the contribution different substrates have on the DBS workflow for the measurement of FAME and the significant contribution of saturated fatty acids in the paper substrate. The synthetic mPPM substrate has a negligible SFA contribution and the mPPM insert for the hemaPEN is ideal for collecting blood volumes from $2.74 - 25 \mu$ L.

Table I. Each class of FAME (n=3) reported as % contribution to the total summed peak area for all FAME (taking into account the relative response of each FAME to the 37 FAME standard). Colored cells represent a difference of >10% from the blood reference.

3 µL	Blood	226 - 6 mm	mPPM
SFA	29.4±0.4%	34.0±0.9%	29.6±0.7%
MUFA	18.3±0.5	18.0±0.4%	19.1±0.2%
 n-6 PUFA	42.5±0.3%	39.0±0.6%	41.7±0.3%
n-3 PUFA	8.3±0.2%	7.6±0.4%	8.3±0.3%
10 µL	Blood	226 - 6 mm	mPPM
SFA	27.7±0.3%	31.9±0.6%	28.0±0.4%
 MUFA	18.2±0.2%	18.4±0.5%	18.2±0.1%
 n-6 PUFA	44.7±0.6%	41.2±0.2%	43.9±0.5%
 n-3 PUFA	8.0±0.1%	6.9±0.1%	8.1±0.1%

Contribution of substrate contamination to % FA measurement



Figure 7. Impact of substrate contamination on the relative abundance of blood sourced fatty acids.

A
FTA Gene
903
GenSaver
GenCollect
FTA Elute
mPPM

B
Image: Study and the state of the state of

Figure 8. 1% Agarose gel electrophoresis of DNA extracts from different DBS paper substrates and mPPM insert. DNA extracts were prepared from A – Qiagen QIAamp DNA and B - Zymo Quick DNA mini.

Conclusion:

All of the DBS substrates resulted in recoveries of DNA although the FTA Elute and mPPM consistently outperformed the other substrates both for yield and high molecular weight quality.





Contact info@trajanscimed.com for further information

The hemaPEN[®] is under development and the prototype device is supplied for research or investigational purposes only. This device is not for therapeutic or diagnostic use. hemaPEN[®] is a registered trademark owned by Trajan Scientific Australia Pty Ltd. PerkinElmer 226™ is a trademark of PerkinElmer, Inc or its affiliated entities. 903[®] is a registered trademark of Eastern Business Forms, Inc. *Claims made in content have not been audited

www.trajanscimed.com