## Mitra® Microsample Collection Kit (RUO)

### IMPORTANT – Read the entire instructions before use!

These instructions apply to the following Mitra® Device configurations:


- Microsample Collection Kit

### Intended Use:

The Microsample Collection Kit is a single-use non-sterile procedure pack intended for the collection, storage and transport of capillary whole blood specimens of specified volume using the contents of the collection kit for Research Use Only (RUO).

### Device Description:

The intended user is a healthcare professional or a layperson (including use with pediatric population). The blood sample is typically taken through a prick in the fingertip using a lancet. The VAMS® sampler tip (contained within the Mitra Device) is gently applied to the drop of blood on the fingertip

until the entire sampler tip turns red indicating a volumetric sample (10, 20, or 30  $\mu$ L depending on the device being used) has been collected. This process is repeated for each sampler tip in the device. Once all sampler tips (one (1), two (2), three (3), or four (4) depending on the device being used) have been filled, the device outer housing is clicked closed to secure the sample(s).

### Materials Provided:

- Mitra Device
- Specimen Bag with Desiccant
- Sample ID Barcode
- Lancet (x2 or x4)
- Gauze (x1 or x2)
- Bandage (x2 or x4)
- Plastic Shipping Envelope

### Product Specifications:

- Sample Type: Dried Whole Blood
- Formats/Configurations: Collection Kit
- Sample Volume: 10, 20, or 30  $\mu$ L
- Substrate: Hydrophilic porous polymer
- Volumetric Precision (%RSD):  $\leq 5\%$
- Typical sampling event (per sampling tip size): 8 sec (10, 20  $\mu$ L) – 12 sec (30  $\mu$ L)
- Number of Samples Collected: Up to 4 (depending on Device used)
- Shelf Life: Refer to expiry date on product label

- Device Storage: Up to 30°C
- Sample Storage: The guidelines for sample storage are analyte dependent and will need to be determined by device end-user

## Warnings and Precautions:

- Single use only for a single individual.
- Do not use after expiration date.
- For external use only.
- Devices should be transported/mailed to the analytical laboratory, and appropriate documentation maintained according to local regulations and the analytical laboratory procedures and policies.
- Do not use if device packaging has been opened or damaged.
- Laboratories must validate use of product for their specific assay.
- Observe universal biological risk precautions.
- All used materials with blood residues must be handled and disposed of safely in accordance with local regulations.
- Any serious incident in relationship with the Mitra Microsample collection kit should be reported as soon as possible to Trajan Scientific and Medical ([neo.support@trajanscimed.com](mailto:neo.support@trajanscimed.com)) and the competent authority of the Member State in which the user and/or patient is established.


**WATCH –**  
Instructional videos at  
[www.neoteryx.com/collect](http://www.neoteryx.com/collect)  
or scan the QR code  
to access.


- Please reach out to [neo.support@trajanscimed.com](mailto:neo.support@trajanscimed.com) in case of multiple under-sampling events

## Instructions for Sampling:

1. Unpack and prepare your kit contents:
  - a. Tear open specimen bag.
  - b. Remove device from specimen bag.  
**IMPORTANT – Leave the desiccant in the specimen bag!**
  - c. Open the device by pulling tabs on device apart from each other until the two flaps are folded downwards. **NOTE:** Do NOT remove sampler tips from sampler bodies. Do NOT remove sampler bodies from device.
  - d. Twist off lancet caps.
  - e. Unwrap gauze.
  - f. Unpeel back of bandages.
2. **Wash hands** with warm, soapy water. Dry hands.
3. **Warm your hands** by rubbing them together. You may also warm your hands by walking, jumping, or any other method you prefer. **Warming your hands is IMPORTANT for good blood flow.**


4. **Select the finger you will use to collect the blood sample.** See illustration for recommended locations off-center of finger midline.

5. **Prick your finger** by laying hand on a hard surface. Position the lancet on the chosen location. Apply firm pressure to end of lancet until a click is heard.
6. **Wipe away first blood drop with the gauze.** A blood drop may take up to 20 seconds to form. If needed, gently massage pricked finger upward until drop forms.

7. **Touch sampler tip to surface of blood sample.**

Watch sampler tip turn FULLY red as it absorbs the blood, count 2 seconds, and remove from blood. It is OK to apply the sampler tip to the blood drop several times to fill until no white is visible.


**IMPORTANT!**

- **Do NOT fully submerge the sampler tip in blood.**
- **Sampler tip should always point downward towards floor as illustrated.**
- **Do NOT drip blood onto the sampler tip.**

8. Repeat step 7 with remaining samplers in device. Device will contain one (1), two (2), three (3), or four (4) samplers depending on device format received.
9. Check the sampler tips (see Figure 1). Upon successful sampling, the entirety of each sampler tip should be colored red with no white visible.

10. **Set device down** on hard surface and apply bandage to your finger.


If you experience excessive bleeding after collection, apply pressure and contact your healthcare provider immediately.

11. **Close device** by lifting the flaps to meet at the top. Press together until a click is heard.
12. **Insert device** into the specimen bag and seal shut. **Ensure the desiccant is still in the bag.**
13. **Insert specimen bag** into plastic shipping envelope. Ensure a pre-addressed label is affixed to mailing envelope. Send to lab for processing as instructed.


14. **Dispose of the lancet** by enclosing it in the plastic pouch that contained the lancets, seal it and trash it in a bin. Do not dispose the lancet in a recycling bin.

## Ensure sampler tips are filled correctly.

Over-sampling occurs when:

1. Blood is dripped onto the sampler tip from above.
2. The entire sampler tip is submerged in a large blood drop formed on the fingertip.

Under-sampling occurs when:

1. The sampler tip is removed from the blood too soon. Touch tip to blood until no white remains.
2. Blood flow stops. If this happens, repeat steps 3-7 with same tip until it turns fully red.
3. Refer to Warnings and Precautions section in case of multiple under-sampling events.

Correctly sampled

Figure 1

**For further assistance, please  
do not hesitate to contact us:  
[neo.support@trajanscimed.com](mailto:neo.support@trajanscimed.com)**

Specifications are subject to change.

Neoteryx®, Mitra® and VAMS® are registered trademarks  
owned by Trajan Scientific Australia Pty Ltd.

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Neoteryx®, Mitra® och VAMS® är registrerade varumärken som ägs av Trajan Scientific Australia Pty Ltd.
**Trajan Scientific Australia Pty Ltd (Sponsor)**
**Axel Semrau GmbH**