Mitra[®] Device (RUO)

EN

DE

DK

ES

FR

IT

NL

NO

PL

PT

Mitra® Device (RUO)

IMPORTANT – Read the entire instructions before use!

These instructions apply to the following Mitra® Device configurations:


- Mitra Device Alone
- Mitra Device in Specimen Bag with Desiccant

Intended Use:

A single-use, non-sterile device used as a specimen collector, and for the storage and transport of blood and other biological fluids for analytical and diagnostic analyses for Research Use Only (RUO).

Device Description:

The intended user is a healthcare professional or a layperson (including use with pediatric population). The blood sample is typically taken through a prick in the fingertip using a lancet. The sampler tip (contained within the Mitra Device) is gently applied to the drop of blood on the fingertip until the entire sampler tip

turns red indicating a volumetric sample (10, 20, or 30 μ L depending on the device being used) has been collected. This process is repeated for each sampler tip in the device. Once all sampler tips (one (1), two (2), three (3), or four (4) depending on the device being used) have been filled, the device outer housing is clicked closed to secure the sample(s).

Materials Provided:

- Mitra Device
- Specimen Bag with Desiccant*
- Sample ID Barcode (if requested)

Materials Required for Blood Sampling, Not Provided by Trajan:

- Disposable Lancet**
- Gauze
- Bandage

* Part #s containing an "A" represent the "alone" configuration and do not contain a specimen bag (e.g. VM-102ANLR). Part #s containing a "S" represent the "specimen bag" configuration and do contain a bag (e.g. VM-102SNLR).

** Use a preferred lancet according to individual's need. Use lancet as per manufacturer's IFU. Contact neo.support@trajanscimed.com if you require assistance in choosing the most appropriate lancet.

Product Specifications:

- Sample Type: Dried Whole Blood and other biological fluids
- Formats/Configurations: Alone or In Specimen Bag with Desiccant
- Sample Volume: 10, 20, or 30 μL
- Substrate: Hydrophilic porous polymer
- Volumetric Precision (%RSD): $\leq 5\%$
- Typical sampling event (per sampling tip size): 8 sec (10, 20 μL) – 12 sec (30 μL)
- Number of Samples Collected: Up to 4 (depending on Device used)
- Shelf Life: Refer to expiry date on product label
- Device Storage: Up to 30°C
- Sample Storage: The guidelines for sample storage are analyte dependent and will need to be determined by device end-user

Warnings and Precautions:

- Single use only for a single individual.
- Do not use after expiration date.
- For external use only.
- Devices should be transported/mailed to the analytical laboratory, and appropriate documentation maintained according to local regulations and the analytical laboratory procedures and policies.
- Do not use if device packaging has been opened or damaged.
- Laboratories must validate use of product for their specific assay.
- Observe universal biological risk precautions.
- All used materials with blood residues must be handled and disposed of safely in accordance with local regulations.
- Please reach out to neo.support@trajanscimed.com for alternative instructions of the collection of other biological fluid. Dried biological fluids can be dried and safely transported as per instructions contained in this IFU.
- Any serious incident in relationship with the Mitra device should be reported as soon as possible to Trajan Scientific and Medical (neo.support@trajanscimed.com) and the competent authority of the Member State in which the user and/or patient is established.
- Please reach out to neo.support@trajanscimed.com in case of multiple under-sampling events.


WATCH –
Instructional videos at
www.neoteryx.com/collect
or scan the QR code
to access.

Recommendations for Finger Blood Sampling


- Warm hands thoroughly (e.g. run hands under hot water, rub hands together, walk briskly for 5 minutes).
- Shake and massage palm and fingers upward to activate blood flow.
- Use middle or ring finger on your non-dominant hand.
- Target placement of the lancet off-center of the finger midline.
- Lay finger you intend to lance on a hard surface while activating the lancet.
- If insufficient blood, use another lancet on a different finger tip.

Instructions for Blood Sampling:


1. **OPTIONAL:** If your device came in a specimen bag, tear open the specimen bag using the pre-slit notch. Remove the device from the specimen bag. Do not remove the desiccant from specimen bag. If your device did not come in a bag, skip to step 2.

2. Open the device by pulling tabs on device apart from each other until the two flaps are folded downwards. **NOTE:** Do NOT remove sampler tips from sampler bodies. Do NOT remove sampler bodies from device.

3. **Touch sampler tip to surface of blood sample (Figure 1).** Watch sampler tip turn **FULLY** red as it absorbs the blood, count 2 seconds, and remove from blood. It is OK to apply the sampler tip to the blood drop several times to fill until no white is visible.


IMPORTANT!

- **Do NOT fully submerge the sampler tip in blood.**
- **Sampler tip should always point downward towards floor as illustrated.**
- **Do NOT drip blood onto the sampler tip.**


4. Repeat step 3 with remaining samplers in device. Device will contain one (1), two (2), three (3), or four (4) samplers depending on device format received.

Ensure sampler tips are filled correctly.

Over-sampling occurs when:

1. Blood is dripped onto the sampler tip from above.
2. The entire sampler tip is submerged in a large blood drop formed on the fingertip.

Under-sampling occurs when:


1. The sampler tip is removed from the blood too soon. Touch tip to blood until no white remains.
2. Blood flow stops. If this happens, repeat steps 3-7 with same tip until it turns fully red.
3. Refer to Warnings and Precautions section in case of multiple under-sampling events.

Correctly sampled

Figure 1

5. Check the sampler tips. Upon successful sampling, the entirety of each sampler tip should be colored red with no white visible.


6. Close device by lifting the flaps to meet at the top. Press together until a click is heard.

7. IF YOU RECEIVED A SPECIMEN BAG: Insert device into the specimen bag and seal shut. Ensure the desiccant is still in the bag.

**For further assistance,
please do not hesitate to contact us:
neo.support@trajanscimed.com**

Specifications are subject to change.

Neoteryx®, Mitra® and VAMS® are registered trademarks owned by Trajan Scientific Australia Pty Ltd.

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Neoteryx®, Mitra® e VAMS® são marcas registadas propriedade da Trajan Scientific Australia Pty Ltd.
Trajan Scientific Australia Pty Ltd (Sponsor)
Axel Semrau GmbH